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A. Absence of gap in the implant-abutment connection

(1) Larrucea Verdugo C, Jaramillo Núñez G, Acevedo Avila A, Larrucea San Martín C. Microleakage of the prosthetic abutment/implant interface with internal and external connection: in vitro study. Clin Oral Implants Res. 2014 Sep;25(9):1078-83.

(2) Larrucea C, Conrado A, Olivares D, Padilla C, Barrera A, Lobos O. Bacterial microleakage at the abutment-implant interface, in vitro study. Clin Implant Dent Relat Res. 2018 Jun;20(3):360-367.

(3) He Y, Fok A, Aparicio C, Teng W. Contact analysis of gap formation at dental implant-abutment interface under oblique loading: A numerical-experimental study. Clin Implant Dent Relat Res. 2019 Aug;21(4):741-752.

(4) Carlos LV, C Carlos N, Sm Karina L, K Sunil B, E Carlos P, G Olga L. Comparative study of bacterial microfiltration in the implant-abutment interface, with straight and conical internal connections, in vitro. Clin Exp Dent Res. 2021 Dec;7(6):1014-1024.

(5) Honório Tonin BS, He Y, Ye N, Chew HP, Fok A. Effects of tightening torque on screw stress and formation of implant-abutment microgaps: A finite element analysis. J Prosthet Dent. 2022 Jun;127(6):882-889.

B. Maintenance of marginal bone

Average bone loss at 5-year follow-up

(6) R. Fernández, D. Peñarrocha-Oltra, M. Peñarrocha-Diogo, J. Buti, E. Khanari, M. Esposito. Natural or platat positioning of immediate post-extractive implants in the aesthetic zone? 5-year results of a multicentre randomised controlled trial. Clinical Trials in Dentistry 2020;02(1):59-71

(7) Serrano B, Sanz-Sánchez I, Serrano K, Montero E, Sanz M. One-year outcomes of dental implants with a hybrid surface macro-design placed in patients with history of periodontitis: A randomized clinical trial. J Clin Periodontol. 2022 Feb;49(2):90-100.

(8) Maceiras L, Liñares A, Nóvoa L, Batalla P, Mareque S, Pérez J, Blanco J. Marginal changes at bone-level implants supporting dental prostheses with or without intermediate standardised abutments after 36 months: Randomised controlled clinical trial. Clin Oral Implants Res.2024; 00:1-9.

(9) Bernabeu-Mira JC, Pellicer-Chover H, Peñarrocha-Diogo M, Peñarrocha-Oltra D. In Vitro Study on Bone Heating during Drilling of the Implant Site: Material, Design and Wear of the Surgical Drill. Materials (Basel). 2020 Apr 19;13(8):1921

C. Survival rate

C.1. Implant survival at 1- and 5-year follow-up

(10) Peñarrocha-Diogo MA, Flichy-Fernández AJ, Ionso-González R, Peñarrocha-Oltra D, Balaguer-Martínez J, Peñarrocha-Diogo M. Influence of implant neck design and implant-abutment connection type on peri-implant health. Radiological study. Clin Oral Implants Res. 2013 Nov;24(11):1192-200.

(11) Pellicer-Chover H, Peñarrocha-Diogo M, Peñarrocha-Oltra D, Gomar-Vercher S, Agustín-Panadero R, Peñarrocha-Diogo M. Impact of crestal and subcrestal implant placement in peri-implant bone: A prospective comparative study. Med Oral Patol Oral Cir Bucal. 2016 Jan 1;21(1):e103-10.

(12) Peñarrocha-Diogo M, Balaguer-Martí JC, Peñarrocha-Oltra D, Balaguer-Martínez JF, Peñarrocha-Diogo M, Agustín-Panadero R. A combined digital and stereophotogrammetric technique for rehabilitation with immediate loading of complete-arch, implant-supported prostheses: A randomized controlled pilot clinical trial. J Prosthet Dent. 2017 Nov;118(5):596-603.

(13) Pellicer-Chover H, Peñarrocha-Oltra D, Aloy-Prosper A, Sanchis-Gonzalez JC, Peñarrocha-Diogo MA, Peñarrocha-Diogo M. Comparison of peri-implant bone loss between conventional drilling with irrigation versus low-speed drilling without irrigation. Med Oral Patol Oral Cir Bucal. 2017 Nov 1;22(6):e730-e736.

(14) Pellicer-Chover H, Peñarrocha-Diogo M, Aloy-Prosper A, Canullo L, Peñarrocha-Diogo M, Peñarrocha-Oltra D. Does Apico-Coronal Implant Position Influence Peri-Implant Marginal Bone Loss? A 36-Month Follow-Up Randomized Clinical Trial. J Oral Maxillofac Surg. 2019 Mar;77(3):515-527.

(15) M. Esposito, L. Sbricoli, J. Buti, U. Uccioli, M. Tallarico. Endodontic retreatment vs dental implants of teeth with an uncertain endodontic prognosis: 5-year results from a randomised controlled trial. Clinical Trials in Dentistry 2020;02(3):27-43.

C.2. Implant survival up to 7-year follow-up in oncologic patients

(16) Cuesta-Gil M, Ochandiano Caicoya S, Riba-García F, Duarte Ruiz B, Navarro Cuéllar C, Navarro Vila C. Oral rehabilitation with osseointegrated implants in oncologic patients. J Oral Maxillofac Surg. 2009 Nov;67(11):2485-96.

(17) Mancha de la Plata M, Gías LN, Díez PM, Muñoz-Guerra M, González-García R, Lee GY, Castrejón-Castrejón S, Rodríguez-Campo FJ. Osseointegrated implant rehabilitation of irradiated oral cancer patients. J Oral Maxillofac Surg. 2012 May;70(5):1052-63.

(18) Navarro Cuellar C, Caicoya SJ, Acero Sanz JJ, Navarro Cuellar I, Muela CM, Navarro Vila C. Mandibular reconstruction with iliac crest free flap, nasolabial flap, and osseointegrated implants. J Oral Maxillofac Surg. 2014 Jun;72(6):1226.e1-15.

(19) Navarro Cuéllar C, Ochandiano Caicoya S, Navarro Cuéllar I, Valladares Pérez S, Fariña Sirandoni R, Antúnez-Conde R, Díez Montiel A, Sánchez Pérez A, López López AM, Navarro Vila C, Salmerón Escobar JJ. Vertical Ridge Augmentation of Fibula Flap in Mandibular Reconstruction: A Comparison between Vertical Distraction, Double-Barrel Flap and Iliac Crest Graft. J Clin Med. 2020 Dec 30;10(1):101.

(20) Antúnez-Conde R, Salmerón JI, Díez-Montiel A, Agea M, Gascón D, Sada Á, Navarro Cuéllar I, Tousidonis M, Ochandiano S, Arenas G, Navarro Cuéllar C. Mandibular Reconstruction With Fibula Flap and Dental Implants Through Virtual Surgical Planning and Three Different Techniques: Double-Barrel Flap, Implant Dynamic Navigation and CAD/CAM Mesh With Iliac Crest Graft. Front Oncol. 2021 Oct 5;11:719712.

(21) Navarro Cuéllar C, Tousidonis Rial M, Antúnez-Conde R, Agea Martínez M, Navarro Cuéllar I, Salmerón Escobar JJ, Navarro Vila C. Functional Outcomes with Facial Artery Musculo-Mucosal (FAMM) Flap and Dental Implants for

Reconstruction of Floor of the Mouth and Tongue Defects in Oncologic Patients. J Clin Med. 2021 Aug 17;10(16):3625.

(22) Navarro Cuéllar C, Tousidonis Rial M, Antúnez-Conde R, Ochandiano Caicoya S, Navarro Cuéllar I, Arenas de Frutos G, Sada Urmeneta Á, García-Hidalgo Alonso MI, Navarro Vila C, Salmerón Escobar JJ. Virtual Surgical Planning, Stereolitographic Models and CAD/CAM Titanium Mesh for Three-Dimensional Reconstruction of Fibula Flap with Iliac Crest Graft and Dental Implants. J Clin Med. 2021 Apr 29;10(9):1922.

D. Bone-to-Implant Contact (BIC) ratio

(23) Blanco J, Alvarez E, Muñoz F, Liñares A, Cantalapiedra A. Influence on early osseointegration of dental implants installed with two different drilling protocols: a histomorphometric study in rabbit. Clin Oral Implants Res. 2011 Jan;22(1):92-9.

(24) Soto-Peñaloza D, Caneva M, Viña-Almunia J, Martín-de-Llano JJ, García-Mira B, Peñarrocha-Oltra D, Botticelli D, Peñarrocha-Diogo M. Effect on osseointegration of two implant macro-designs:A histomorphometric analysis of bicortically installed implants in different topographic sites of rabbit's tibiae. Med Oral Patol Oral Cir Bucal. 2019 Jul 1;24(4):e502-e510.

(25) Monje A, Chappuis V, Monje F, Muñoz F, Wang HL, Urban IA, Buser D. The Critical Peri-implant Buccal Bone Wall Thickness Revisited: An Experimental Study in the Beagle Dog. Int J Oral Maxillofac Implants. 2019 November/December;34(6):1328-1336.

(26) Sanchez-Perez A, Cachazo-Jiménez C, Sánchez-Matás C, Martín-de-Llano JJ, Davis S, Carda-Batalla C. Effects of Ultraviolet Photoactivation on Osseointegration of Commercial Pure Titanium Dental Implant After 8 Weeks in a Rabbit Model. J Oral Implantol. 2020 Apr 1;46(2):101-107.

E. Implant stability Values (ISQ)

(27) Gómez-Polo M, Ortega R, Gómez-Polo C, Martín C, Celemin A, Del Río J. Does Length, Diameter, or Bone Quality Affect Primary and Secondary Stability in Self-Tapping Dental Implants? J Oral Maxillofac Surg. 2016 Jul;74(7):1344-53.

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